

Ionode Pty Ltd  
**SAFETY DATA SHEET**

Date Prepared: November, 2019  
Version No: 1.2

---

**1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER**

---

Product Name: Buffer Solution pH10 Colour Coded Blue  
Product Codes: pH10.00-250 and pH10.00-1L  
Other Names: Nil  
Uses: Analytical Reagent

Supplier: Ionode Pty Ltd  
12 Walker Street, Tennyson Qld 4105

Contacts: Telephone: 61 07 38481660  
Emergency Phone: 61 07 38481660

---

**2. HAZARDS INFORMATION**

---

**Hazard classification:** Non Hazardous. Non Dangerous Goods.

**Risk phrases:**

Not considered a hazard according to the criteria of NOHSC.

**Safety phrases:**

Not considered a hazard according to the criteria of NOHSC.

---

**3. COMPOSITION / INFORMATION ON INGREDIENTS**

---

**Ingredients :**

Chemical Entity	CAS No	Proportion
Sodium tetraborate	[1330-43-4]	<10%
Sodium hydroxide	[1310-73-2]	<1%
Water	[7732-18-5]	to 100%

---

**4. FIRST AID MEASURES**

---

Safety showers and eye wash facilities should be provided.

**Swallowed :**

If conscious wash out mouth with water. Seek medical advice. Show this MSDS to medical practitioner.

**Eye :**

Immediately hold eyelids open and flood with water for at least 15 minutes. Obtain medical aid. Show this MSDS to medical practitioner.

**Skin :**

Remove contaminated clothing. Immediately wash skin thoroughly with water and mild soap. Seek medical advice if irritation persists. Show this MSDS to medical practitioner. Launder clothing before reuse.

**Inhaled :**

Remove from contaminated air. Maintain breathing with artificial respiration if necessary. Seek medical assistance. Show this MSDS to a doctor.

---

## 5. FIRE FIGHTING MEASURES

---

### Suitable Extinguishing Media:

Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.

### Hazards from Combustion Products:

Product will not burn or support combustion. Decomposition products include oxides of sodium and boron.

### Precautions for Fire Fighters and Special Protective Equipment:

Fire fighters and others who may be exposed to combustion products during fire should wear full protective clothing including positive pressure self-contained breathing apparatus (SCBA). Wear SCBA with full face-piece, operated in positive pressure mode when fighting fires.

---

## 6. ACCIDENTAL RELEASE MEASURES

---

### Emergency procedures:

Prevent from entering waterways. Restrict access to area. Ventilate area. Remove chemicals that can react with the spilled material.

### Methods and materials for containment and clean up:

Use inert material such as sand or earth to contain spill or leak. Absorb spills with chemical absorber or vermiculite and dispose of in accordance with local regulations.

---

## 7. HANDLING AND STORAGE

---

### Precautions for Safe Handling:

Do not get in eyes, on skin, on clothing. Avoid prolonged or repeated exposure.

### Conditions for Safe Storage:

Store sealed in original container in a cool well ventilated situation away from foods and other chemicals. Do not store in direct sunlight. Observe good hygiene and housekeeping practices.

---

## 8. EXPOSURE CONTROLS / PERSONAL PROTECTION

---

### National Exposure Standards:

SWA Australia – Borates, tetra, sodium salts (decahydrate) – 5mg/m<sup>3</sup>  
Sodium hydroxide – 2mg/m<sup>3</sup> & Peak limitation

**Biological Limit Values:** No data available.

### Engineering Controls:

Not required with normal use.

### Personal Protective Equipment (PPE):

The use of nitrile or neoprene gloves complying with AS 2161 and the use of faceshield, chemical goggles or safety glasses with side shield protection complying with AS/NZS 1337 is recommended.

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

---

<b>Appearance :</b>	Clear blue liquid
<b>Odour:</b>	Nil
<b>pH:</b>	10
<b>Boiling Point (°C) :</b>	Not applicable
<b>Freezing/melting Point:</b>	Not applicable
<b>Vapour Pressure (mm of Hg @ 25°C) :</b>	Not applicable
<b>Vapour Density:</b>	Not applicable
<b>Specific Gravity :</b>	1
<b>Flash Point (°C) :</b>	Not flammable
<b>Flammability Limits (%) :</b>	Not flammable
<b>Solubility in Water (g/L) :</b>	Soluble

---

## 10. STABILITY AND REACTIVITY

---

### Chemical stability:

Stable.

### Conditions to avoid:

Excessive heat. Strong Sunlight, Absorption of carbon dioxide

### Incompatible materials:

Acids, alkalis

### Hazardous decomposition products:

Refer to section 5 (Fire Fighting Measures).

### Hazardous reactions:

Hazardous polymerization will not occur.

---

## 11. TOXICOLOGICAL INFORMATION

---

### Health Effects:

**Swallowed** : May cause irritation of the gastric system. Ingestion of large quantities may cause severe vomiting, diarrhea, shock or death. For sodium tetraborate LD<sub>50</sub> : oral infant 1000mg/kg, oral man 709mg/kg LD50 oral rat 2660mg/kg

**Eye** : May be irritating to eye tissue. For sodium hydroxide 500mg applied to rabbit skin produced severe irritation after 24 hours

**Skin** : May be irritating to skin tissue. May be harmful by skin absorption.

**Inhaled** : Not considered a hazard with normal laboratory use. Mists may cause irritation of mucous membranes.

**Chronic Effects**: No data available.

---

## 12. ECOLOGICAL INFORMATION

---

### Ecotoxicity:

No data available.

### Persistence and degradability:

No data available.

### Mobility:

No data available.

---

## 13. DISPOSAL CONSIDERATIONS

---

Contact a licensed professional waste disposal service to dispose of this material. Observe all federal, state and local environmental regulations.

---

## 14. TRANSPORT INFORMATION

---

**UN Number:** None allocated

**UN Proper Shipping Name:** None allocated

**Class and subsidiary risk(s):** None allocated

**Packing Group:** None allocated

**Hazchem Code:** None allocated

**Special precautions for user** : Nil

---

## **15. REGULATORY INFORMATION**

---

**Standard for the Uniform Scheduling of Drugs and Poisons (SUSDP):**

Not scheduled

---

## **16. OTHER INFORMATION**

---

**Disclaimer:**

All information given by the Company is offered in good faith and is believed to the best of our knowledge to be accurate. However this information is offered without warranty representation inducement or licence and the Company does not assume legal responsibility for reliance upon the same.

Every person dealing with the materials referred to herein does so at his or her own risk absolutely and must make independent determinations of suitability and completeness of information from all sources to ensure their proper use.