



IJ40A-HT - Intermediate Junction High Temperature pH Electrode



SCALE 1.5:1

Serviceable 'Intermediate Junction'
(Secondary Junction) *



Slim Junction Sleeve Provided

Temp. Range: 0-100°C

Body Material: Kynar (PVDF)

pH Range: 0-14

Reference Potential: -44mV vs SCE

Junction: Sleeve Junction

Max. Pressure Rating: 40 PSI

Submersible: Yes

ionode

GOOD CHEMISTRY



*3M KCl Gel Electrolyte supplied. User may fill with their application specific electrolyte.