

Ionode Pty Ltd
SAFETY DATA SHEET

Date Prepared: November, 2019
Version No: 1.1

1. IDENTIFICATION OF THE MATERIAL AND SUPPLIER

Product Name: Potassium Nitrate Reference Gel Electrolyte
Product Codes: KN45G
Other Names: Nil
Uses: Analytical Reagent

Supplier: Ionode Pty Ltd
12 Walker Street, Tennyson Qld 4105

Contacts: Telephone: 61 07 38481660
Emergency Phone: 61 07 38481660

2. HAZARDS INFORMATION

Hazard classification: Classified as non-hazardous according to the Globally Harmonised System of classification and labelling of Chemicals (GHS) including Work, Health and Safety regulations, Australia

3. COMPOSITION / INFORMATION ON INGREDIENTS

Ingredients :

Chemical Entity	CAS No	Proportion
Potassium nitrate	[7757-79-1]	10%
Triethylene glycol	[112-27-6]	>10%
Hydroxyethyl cellulose	[9004-62-0]	>2%
Water	[7732-18-5]	to 100%

4. FIRST AID MEASURES

Safety showers and eye wash facilities should be provided.

Swallowed :

If conscious wash out mouth with water. Seek medical advice. Show this MSDS to medical practitioner.

Eye :

Immediately hold eyelids open and flood with water for at least 15 minutes. Obtain medical aid. Show this MSDS to medical practitioner.

Skin :

Remove contaminated clothing. Immediately wash skin thoroughly with water and mild soap. Seek medical advice if irritation persists. Show this MSDS to medical practitioner. Launder clothing before reuse.

Inhaled :

Remove from contaminated air. Maintain breathing with artificial respiration if necessary. Seek medical assistance. Show this MSDS to a doctor.

5. FIRE FIGHTING MEASURES

Suitable Extinguishing Media:

Water spray. Carbon dioxide, dry chemical powder, or appropriate foam.

Hazards from Combustion Products:

Solution will not burn or support combustion. Decomposition products include oxides of carbon and potassium.

Precautions for Fire Fighters and Special Protective Equipment:

Fire fighters and others who may be exposed to combustion products during fire should wear full protective clothing including positive pressure self-contained breathing apparatus (SCBA).

6. ACCIDENTAL RELEASE MEASURES

Emergency procedures:

Prevent from entering waterways. Restrict access to area. Remove chemicals that can react with the spilled material.

Methods and materials for containment and clean up:

Wear protective clothing when dealing with spills. Use inert material such as sand or earth to contain spill or leak. Absorb spills with chemical absorber or vermiculite and dispose of in accordance with local regulations.

7. HANDLING AND STORAGE

Precautions for Safe Handling:

Do not get in eyes, on skin, on clothing. Avoid prolonged or repeated exposure.

Conditions for Safe Storage:

Store sealed in original container in a cool well ventilated situation away from foods and other chemicals. Do not store in direct sunlight. Observe good hygiene and housekeeping practices.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

National Exposure Standards:

Worksafe – None known

Biological Limit Values: No data available.

Engineering Controls:

Not required with normal use. If mists are likely to be generated maintain atmospheric concentrations to minimal levels with extraction ventilation.

Personal Protective Equipment (PPE):

The use of nitrile or neoprene gloves complying with AS 2161 and the use of faceshield, chemical goggles or safety glasses with side shield protection complying with AS/NZS 1337 is recommended.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance :	Clear gel
Odour:	Nil
pH:	7
Boiling Point (°C) :	Not known
Freezing/melting Point:	Not known
Vapour Pressure (mm of Hg @ 25°C) :	Not known
Vapour Density:	Not applicable
Specific Gravity :	Not known
Flash Point (°C) :	Not flammable
Flammability Limits (%) :	Not flammable
Solubility in Water (g/L) :	Soluble

10. STABILITY AND REACTIVITY

Chemical stability:

Stable.

Conditions to avoid:

Excessive heat.

Incompatible materials:

Materials that react with water

Hazardous decomposition products:

Refer to section 5 (Fire Fighting Measures).

Hazardous reactions:

Hazardous polymerization will not occur.

11. TOXICOLOGICAL INFORMATION

Health Effects:

Swallowed : May be irritating to tissue. May be harmful by if swallowed. Ingestion of large quantities may cause vomiting, diarrhoea. Ingestion of large quantities may cause heart condition due to high potassium level. For triethylene glycol LDLo Oral Human 5000 mg/kg

Eye : May be irritating to eye tissue.

Skin : May be irritating to skin tissue. May be harmful by skin absorption.

Inhaled : May be irritating to respiratory tissue.

Chronic Effects: May be harmful with repeated or prolonged exposure.

12. ECOLOGICAL INFORMATION

Ecotoxicity:

No data available.

Persistence and degradability:

No data available.

Mobility:

No data available.

13. DISPOSAL CONSIDERATIONS

Contact a licensed professional waste disposal service to dispose of this material. Observe all federal, state and local environmental regulations.

14. TRANSPORT INFORMATION

UN Number: Not applicable

UN Proper Shipping Name: Not applicable

Class and subsidiary risk(s): Not applicable

Packing Group: Not applicable

Hazchem Code: Not applicable

Special precautions for user : Nil

15. REGULATORY INFORMATION

Standard for the Uniform Scheduling of Drugs and Poisons (SUSDP):

Not Scheduled

16. OTHER INFORMATION

Disclaimer:

All information given by the Company is offered in good faith and is believed to the best of our knowledge to be accurate. However this information is offered without warranty representation inducement or licence and the Company does not assume legal responsibility for reliance upon the same.

Every person dealing with the materials referred to herein does so at his or her own risk absolutely and must make independent determinations of suitability and completeness of information from all sources to ensure their proper use.